

D-B-Life

Infant reusable resuscitator



The image can differ from the actual product



Reusable self inflating manual resuscitator set for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

Connection standard ISO 15/22 mm
Possibility to connect the bag to sources of oxygen supply.
Steam sterilizable
Suitable for patients up to 10 kg approx.

Included equipment

Reservoir
Overp. and non-rebreathing valve
Infant mask (sz1)
Anti crush oxygen tube

Technical data

Materials	Silicone, PVC, PC
Bag volume	340 ml
Reservoir volume	1600 ml
Stroke volume range	150 ml
Overpressure valve	40 cmH ₂ O
Dead space	7ml
Inspirator and expiratory resistance at 50l/min	< 5 cmH ₂ O
Overall weight	350 g

Class IIa med. Dev. according to Dir. 93/42/CEE

D - B-LIFE INFANT MANUAL
RESUSCITATOR SET
BG03202



Distributed by Spencer Italia s.r.l.

Rev.0 (29/03/2021)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it