

B-BAG

Double bag resuscitator



Double bag self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The double bag transmits the pressure variations to the operator more effectively
 Possibility to connect the bag to sources of oxygen supply
 Low storage space required
 Latex Free mask
 Steam sterilizable

Technical data

| | |
|--------------------|---------------------------|
| Materials | Rubber, Silicone, PC, PVC |
| Bag volume | 1600 ml |
| Reservoir volume | 2600 ml |
| Overpressure valve | 60 cmH ₂ O |
| Weight of the kit | 570 g |

Standard equipment

| | |
|---------|---------------------------------|
| BG07002 | Res.bag w/valve for B-Bag adult |
| EV50011 | EVX 11 – Non reb. Overpr. valve |
| RM20705 | Silicone face mask sz.5 |

Med. Dev. CE marked according to Dir. 93/42/CEE

B-BAG - DOUBLE BAG W/OVERP.
 RES. SP. MASK

BG07000B

Registration number



58824

Optionals

| | |
|----------|---------------------------------------|
| OX20131A | OXX30 anti crushing oxygen tube 10pcs |
|----------|---------------------------------------|

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it