

## B-BAG

### Double bag resuscitator



Double bag self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

#### Specific features

The double bag transmits the pressure variations to the operator more effectively  
 Possibility to connect the bag to sources of oxygen supply  
 Low storage space required  
 Latex Free mask  
 Steam sterilizable

#### Technical data

Materials	Rubber, Silicone, PC, PVC
Bag volume	1600 ml
Reservoir volume	2600 ml
Overpressure valve	60 cmH <sub>2</sub> O
Weight of the kit	570 g

#### Standard equipment

BG07002	Res.bag w/valve for B-Bag adult
EV50011	EVX 11 – Non reb. Overpr. valve
RM20705	Silicone face mask sz.5

Med. Dev. CE marked according to Dir. 93/42/CEE

B-BAG - DOUBLE BAG W/OVERP.  
 RES. SP. MASK

**BG07000B**

Registration number



58824

#### Optionals

OX20131A	OXX30 anti crushing oxygen tube 10pcs
----------	---------------------------------------

Rev.1 (03/04/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)