

CO B-Life PVC Resuscitator



Disposable self inflating manual resuscitator kit for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
Intended use ADULT marked on the device to limit risk of incorrect applications
Connection standard ISO 15/22 mm
Low storage space required
Latex Free

Technical data

Materials	Medical grade PVC, PC
Bag Volume	1600 ml
Reservoir volume	2500 ml
Overpressure valve	60 cmH ₂ O
Maximum oxygen concentration	95%
Max. re expansion rhythm	45 bpm

Standard equipment

BG03091A	Sacca Reservoir
EV50009C	EVX09 Overpressure/non-rebreathing valve
RM10892A	Cuffed mask Size 3
RM10893	Cuffed mask Size 4
OX20130A	OXX30 anti crush oxygen tube

MD CE marked according to Dir. 93/42/CEE

CO B-LIFE - ADULT
W/OVERP.RES.MASK.T.OSS. R-
VALVE

BG08004A



Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it