

## CO B-Life PVC Resuscitator



Disposable self inflating manual resuscitator kit for mechanical ventilation of trauma patients or with breathing difficulties.

### Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands  
 Intended use INFANT marked on the device to limit risk of incorrect applications  
 Connection standard ISO 15/22 mm  
 Low storage space required  
 Latex Free

### Technical data

Materials	Medical grade PVC, PC
Bag Volume	280 ml
Overpressure valve	40 cmH <sub>2</sub> O
Maximum oxygen concentration	95%
Max. re expansion rhythm	95 bpm

### Standard equipment

EV50008C	EVX08 Overpressure/non-rebreathing valve
----------	--

MD CE marked according to Dir. 93/42/CEE

CO B-LIFE - PVC INFANT W/VALVE  
 OVERP. AND R-VALVE  
 BG08010A



### Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)