

CO B-Life PVC Resuscitator



Disposable self inflating manual resuscitator kit for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
Intended use PEDIATRIC marked on the device to limit risk of incorrect applications
Connection standard ISO 15/22 mm
Low storage space required
Latex Free

Standard equipment

BG03091A	Reservoir
EV50008C	EVX08 Overpressure/non-rebreathing valve
RM10890A	Cuffed mask Size 1
RM10891A	Cuffed mask Size 2
OX20130A	OXX30 anti crush oxygen tube
EV50110K	Disposable reanimation filter

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Technical data

Materials	Medical grade PVC, PC
Bag Volume	500 ml
Reservoir volume	2500 ml
Overpressure valve	40 cmH ₂ O
Maximum oxygen concentration	95%
Tube length	2 m
Max. re expansion rhythm	100 bpm

MD CE marked according to Dir. 93/42/CEE

CO B-LIFE -
PA.RES.W/OVERP.RES.MASK.T.OX.
R-V FILT

BG08012A



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it