

CO B-Life PVC Resuscitator



Disposable self inflating manual resuscitator kit for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Intended use INFANT marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Low storage space required
 Latex Free

Technical data

Materials	Medical grade PVC, PC
Bag Volume	280 ml
Reservoir volume	600 ml
Overpressure valve	40 cmH ₂ O
Maximum oxygen concentration ^(with reservoir)	95%
Max. re expansion rhythm	95 bpm

Standard equipment

RM10889A	CUFFED MASK SZ.0 W/OUT VALVE
RM10890A	CUFFED MASK SZ.1 W/OUT VALVE
EV50008C	EVX08 Overpressure/non-rebreathing valve
OX20130A	OXX30 – Anti-crush oxygen tube
BG03092A	Reservoir

MD CE marked according to Dir. 93/42/CEE

CO B-LIFE - INFANT RESUSC.
 W/OVERPR. VALVE W/RES.
 CUFFED MASKS OXYGEN TUBE
 AND R-VALVE
BG08107A



Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it