

ALTO MAX TF

Stretcher with hydraulic system



Spencer ALTO MAX is a sturdy stretcher designed for transport of patients inside hospitals. It is characterized by a hydraulically operated height adjustment system that confers high versatility to the device and minimizes fatigue.

Specific Features

Polyethylene board easy to clean
Adjustable backrest
Trendelenburg and Fowler antishock positions
4 pivoting wheels with brake 2 of which are semi-locked
Foldable side bars <i>Safe bar</i>
Equipped with basket for storage of tools or devices
Oxygen Tank carrier installed on the main frame, suitable 5/7/10 liters oxygen tanks
Additional anterior oxygen tank carrier
Foot operated hydraulic height adjustment system
Can be shortened for a better handling in confined spaces
Integrated foldable IV holder

Standard Equipment

STX 702 two pieces black reflex metal belts

Technical data

Length	2000 ± 10mm
Length folded	1640 ± 20mm
Width	680 ± 10mm
Weight	73 ± 2 kg
Loading capacity	250 kg
Wheels	ø200 mm
Patient floor height	From 650 to 1060 ±20 mm
Basket dimensions	450x640x80 ± 10 mm
Tank container dim.	ø145x740 ± 10 mm (internal)
Anterior tank carrier dim.	ø110x340 ± 20mm (internal)
Frame material	Steel
Side bars length	680 mm
Side bars height	200 mm

Accessories

ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Class I Med.Device according to Dir. 93/42/CEE

ALTO MAX TF YELLOW BOARD
CHROME FRAME W/ HYDRAULIC
SYSTEM



BO00030C

CND Classification	V9099
Registration number	89535

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it