

## BOB

### Stretcher support



Stretcher support intended for installation on ambulances. The device can be extracted to facilitate both loading and unloading of the stretcher and, thanks to the sideways movement, allow easy access to the stretcher.

#### Specific features

- The patented design increases the available space for seat and walk of rescuers
- The structure of the support plane makes it non-deformable
- Loading system with roller (slide optional)
- The control panel can be integrated in the structure and its position can be customized
- Can assume tilt front and tilt back positions
- Dedicated compartment to store immobilization systems like spinal boards or pick-up stretchers.
- Aluminum platform easy to clean
- Electrically operated tilt front and tilt back positions
- Sideways movable and extractable
- Operating handles are easily accessible

#### Accessories

- BO10103 – Electric pack (Led lights and 12V cig. lighter socket)
- BO10100 – Slide 120 mm
- BO10101 – Slide 180 mm
- BO10102 – Slide 350 mm
- Multi-purpose scissors holder integrated in the slide
- Oxygen socket on the front side
- Side Oxygen module (up to 6 sockets)

#### Technical data

Length <small>(Control panel not considered)</small>	2010 mm
Width	600 mm
Platform height	280 mm
Sideways movement (9 steps)	230 ± 10 mm
Materials	Al, Steel, PE, Nylon
Loading capacity	300 kg
Weight	75 ± 10 kg

#### Reference standard

10g tested according to EN1789:2007+A2:2014

#### Device in compliance with dir. 2006/42/CE

BO10005 - BOB – ELECTR. STRETCHER  
SUPP. EXTR. W/SIDEWAYS MOVEM.



Rev.1 (15/06/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)