


BOB 00

Stretcher support



Stretcher support intended for installation on ambulances.

Specific features

Characterized by low weight and essential structure, allows an easy installation
 The structure of the support plane makes it non-deformable
 Loading system with roller (slide optional)
 Extruded aluminum platform easy to clean
 The device allows to increase the rigidity of the stretcher area and makes the loading of the stretcher easier if a proper slide is installed
 Supplied with 180mm ramp

Reference standard

10g tested according to EN1789:2007+A2:2014

Technical data

Length	1950 mm
Width	610 mm
Platform height	55 mm
Materials	Al,
Loading capacity	300 kg
Weight	25 ± 2 kg

BO10010 - BOB 00 - ULTRA THIN
 STRETCHER SUPPORT NO HOUSING

Rev.1 (13/11/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it