


## BOB SLIM

### Stretcher support



Stretcher support intended for installation on ambulances.

#### Specific features

Characterized by low weight and essential structure, allows an easy installation  
 The structure of the support plane makes it non-deformable  
 Loading system with roller (slide optional)  
 Extruded aluminum platform easy to clean  
 The device allows to increase the rigidity of the stretcher area and makes the loading of the stretcher easier if a proper slide is installed

#### Reference standard

10g tested according to EN1789:2007+A2:2014

#### Accessories

BO10103 - Electric pack (Led light and 12V cig. lighter socket)  
 BO10100 – Slide 120 mm  
 BO10101 – Slide 180 mm  
 BO10102 – Slide 350 mm  
 Multi-purpose scissors holder integrated in the slide  
 Oxygen socket on the front side  
 Side Oxygen module (up to 6 sockets)

#### Technical data

Length	1950 mm
Width	610 mm
Platform height	55 mm
Materials	Al,
Loading capacity	300 kg
Weight	25 ± 2 kg

BO10010 - BOB SLIM - ULTRA THIN  
 STRETCHER SUPPORT NO HOUSING

Rev.0 (02/08/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)