


BOB

Electric pack



Optional led lights and 12V cigarette lighter socket to be factory installed on BOB stretcher supports.

Specific features

- The light can be turned on or off if required
- 12V cigarette lighter socket with cover

Technical data

Led bar length	1825 mm
Cigarette lighter output	12V – max 20A

BO10103 - BOB - ELECTRIC PACK

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it