


BOB

Ramp 500 mm



Optional slide to be installed on BOB stretcher supports to allow an easier loading of the stretcher.

Specific features

- Stainless steel body to increase resistance
- Can be raised to allow access to the stretcher support housing
- The stretcher support shall be equipped with dedicated adjustment pads

Technical data

Length	500 mm
Material	Stainless steel

BO10104 - 500mm RAMP FOR BOB

Rev.0 (06/12/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it