


BOB

External control panel



This accessory allows to free the surface of the stretcher support. It can be installed where preferred inside the vehicle according to the length of the connection cable.

Specific features

Easy to install

Has same features as the original display installed on the stretcher support

Technical data

Cable length	5 m
Control panel dimensions	Approx 150x75x35 mm

BO10105 - BOB EXTERNAL CONTROL
PANEL FOR GALAXY

Rev.0 (12/06/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it