

BOB

Stretcher support



Stretcher support with double housing intended for installation on ambulances.

Specific features

The patented design increases the available space for seat and walk of rescuers
 The structure of the support plane makes it non-deformable
 Loading system with roller (slide 120/180/300mm optional)
 Two dedicated compartments to store immobilization systems like spinal boards or pick-up stretchers.
 Extruded aluminum platform easy to clean
 Tested 10g in accordance to EN1789:2007+A2:2014 and compatible with fasteners in compliance with EN 1865 and EN 1789

Accessories

BO10103-Electric pack (Led light and 12V cigarette lighter socket)
 Multi-purpose scissors holder integrated in the slide
 Patented system for equipment with suction devices, oxygen administration and respiratory equipment integrated in the main structure, allowing patient support and increasing convenience for operators
 Oxygen socket on the front side
 Side Oxygen module (up to 6 sockets)

Technical data

Length	2050 mm
Width	600 mm
Platform height	210 mm
Materials	Al, Steel, PE, Nylon
Loading capacity	300 kg
Weight	42 ± 10 kg

BO10201 - BOB - FIXED STRETCHER
 SUPPORT W/DOUBLE HOUSING

Rev.1 (18/07/17)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it