

# BOB 20G

## Stretcher support



Stretcher support intended for installation on ambulances.

### Specific features

The patented design increases the available space for seat and walk of rescuers  
 The structure of the support plane makes it non-deformable  
 Loading system with roller (slide optional)  
 Dedicated compartment to store immobilization systems like spinal boards or pick-up stretchers.  
 Extruded aluminum platform easy to clean

### Technical data

Length	2010 mm
Width	600 mm
Platform height	200 mm
Materials	Al, Steel, PE, Nylon
Loading capacity	300 kg
Weight	42 ± 10 kg

### Reference standard

20g tested according to EN1789:2007+A2:2014

### Accessories

BO10103 - Electric pack (Led light and 12V cig. lighter socket)  
 BO10100 – Slide 120 mm  
 BO10101 – Slide 180 mm  
 BO10102 – Slide 350 mm  
 Multi-purpose scissors holder integrated in the slide  
 Oxygen socket on the front side  
 Side Oxygen module (up to 6 sockets)

BO20001 – BOB 20G - FIXED  
 STRETCHER SUPPORT

Rev.0 (15/06/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)