

BOB 20G

Stretcher support



Stretcher support intended for installation on ambulances. The device can be extracted to facilitate both loading and unloading of the stretcher and, thanks to the sideways movement, allow easy access to the stretcher.

Specific features

- The patented design increases the available space for seat and walk of rescuers
- The structure of the support plane makes it non-deformable
- Loading system with roller (slide optional)
- Dedicated compartment to store immobilization systems like spinal boards or pick-up stretchers.
- Extruded aluminum board easy to clean
- Operating handles are easily accessible

Technical data

Length	2010 mm
Width	600 mm
Support surface height	210 mm
Extraction	670 ± 30 mm
Sideways movement	300 ± 10 mm
Materials	Al, Steel, PE, Nylon
Loading capacity	300 kg
Weight	55 ± 10 kg

Accessori

- BO10103 – Electric pack (Led lights and 12V cig. lighter socket)
- BO10100 – Slide 120 mm
- BO10101 – Slide 180 mm
- BO10102 – Slide 350 mm
- Multi-purpose scissors holder integrated in the slide
- Front oxygen socket
- Side oxygen distribution system (up to 6 sockets)

Reference standard

20g tested according to EN1789:2007+A2:2014

Device in compliance with dir. 2006/42/CE

BO20004 - BOB 20G MECHANICAL
STRETCHER SUPPORT WITH SIDEWAYS
MOVEMENT AND EXTRACTABLE



Rev.0 (15/06/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it