

Burn B BS

Burn Gel



Gel for first aid emergency burns. Used by ambulance services, hospital emergency departments, medical clinics and workplaces and supplied in first aid kits as a first-response treatment for burns.

Specific features

Helps to protect the wound from contaminations, keeping it moist and clean.

Cools the burn and reduce pain

Absorbs and dissipates the heat

Non-toxic

Always ready for use

Shelf life 5 years from manufacturing date

Technical data

Dimensions	85x25x170 mm
Gel	Water based solution

Class IIb Med.Dev. according to Dir. 93/42/CEE

BURN BLOT GEL PACK ML

3.5

BS04000D

CND Classification



M04040602

Distributed by Spencer Italia s.r.l.

Rev.0 (18/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it