

Burn S

Burn stop dressing



These sterile gel- impregnated first aid emergency burn dressings can be used for all burns. They are used by ambulance services, hospital emergency departments, medical clinics and workplaces and supplied in first aid kits as a first-response treatment for burns.

Specific Characteristics

The sterile dress helps to protect the wound and prevent contamination
 Exerts a cooling and soothing action
 Non-adherent to allow removal without pain
 Non-toxic
 Always ready for use

Technical data

Dimensions	10x10 cm
Gel composition (min 35g)	90% water 1% Melaleuca Oil Surfactant and emulsifier
Material	Open cell foam

Class IIa Med.Device according to Dir. 93/42/CEE

BURN S (BA) 10 X 10 CM -
 BURN STOP DRESSING

BS04011D

CND Classification



M04040502

Registration number

74356

Distributed by Spencer Italia S.r.l.

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it