

# Burn S

## Burn stop dressing



These sterile gel- impregnated first aid emergency burn dressings can be used for all burns. They are used by ambulance services, hospital emergency departments, medical clinics and workplaces and supplied in first aid kits as a first-response treatment for burns.

### Specific features

- The sterile dressing helps to protect the wound from contaminations, keeping it moist and clean.
- Cools the burn and reduce pain
- Absorbs and dissipates the heat
- Non-adherent to allow painless removal
- Non-toxic
- Always ready for use
- Shelf life 5 years from manufacturing date

### Technical data

Dimensions	20x20 cm
Gel	Water based solution
Material	Open cell foam

Class IIb Med.Dev. according to Dir. 93/42/CEE

BURN S BA 20 - BURN STOP  
20X20 cm

BS04020D



CND Classification M04040602

Registration number 55180

Distributed by Spencer Italia s.r.l.

Rev.0 (03/10/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)