

BURN C

Fire blanket



This blanket can be used as protection for rescue workers entering or leaving a fire area. It can also be used to extinguish flames on a burn victim and, if necessary, be used as a large dressing, providing instant pain relief and beginning the healing process while victims are being transported to emergency care units.

Specific Characteristics

- Sterile trauma blanket
- Shoulder strap make it easy to carry
- Suitable for: extinguish flame, large dressing, entry or exit from fire
- No loose fibers
- Relieves pain
- Cools and provides the necessary water supply to the skin
- Non-adherent and non-toxic

Technical data

| | |
|---------------------|--|
| Dimensions | 160x120 cm |
| Idrogel composition | 96% water, Melaleuca Alternifolia, emulsifiers |
| PH Gel | 5,5-7 |
| Sterilized by | Gamma ray |
| Blanket Material | Wool |

Class IIb MD according to Dir. 93/42/CEE

BURN C 160 X 120 CM - FIRE BLANKET
BS04050



| | |
|---------------------|---------|
| CND Classification | M040499 |
| Registration number | 67349 |

Distributed by Spencer Italia S.r.l.

Rev.0 (27/09/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it