



BOB MODELS







MECHANICAL STRETCHER SUPPORTS											


Model	BOB SLIM (BOXXXXX)	BO10001	BO10002	BO10003	BO10004	BO10006 – COMPAS M
General specs	Fixed stretcher support	Fixed stretcher support	Sideways movable	Extractable	Sideways movable and extractable	Sideways movable and extractable With tilt front and tilt back position
Surface height	5,5 cm	20 cm	21 cm	21 cm	21 cm	37,5 cm
Loading point when extracted	-	-	-	-	-	From 13 to 37,5 cm above installation floor
Sideways movement	-	No	300 mm	No	300 mm	230 mm
Extraction tilted	-	No	No	No	No	Yes
Housing	-	Yes	Yes	Yes	Yes	Yes
Double housing version 	NA	BO10201	BO10202	BO10203	BO10204	BO10206 (CODE TO BE CREATED)

OPTIONAL						
Slide 	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm
Oxygen pack 	NA	CODE TO BE CREATED	2 Sockets BS – OX07898	CODE TO BE CREATED	CODE TO BE CREATED	CODE TO BE CREATED
Electric pack 	NA	BO10103	BO10103	BO10103	BO10103	BO10103









Model	GALAXY - (BO10000)			MASTER - (BO10008)			LINK - (BO10009)		
General specs	Sideways movable and extractable With tilt front and tilt back position			Sideways movable and extractable			Sideways movable and extractable With tilt back position		
Surface height	Adjustable from 37,5 to 58cm			Adjustable from 28 to 37 cm			Adjustable from 25 to 34cm		
Loading point when extracted	From 16 to 37,5 cm above installation floor			From 10cm below to 28 cm above installation floor			From 4,5 cm below to 25 cm above installation floor		
Sideways movement	230 mm			230 mm (no lateral operation lever)			230 mm (no lateral operation lever)		
Forward movement	No			Yes (150 mm)			No		
Extraction tilted	Yes			Yes			Yes		
Suspension	Yes			Yes			Yes		
Housing	Yes			Yes			No		
Double housing version	 BO10200			BO10208 (TO BE CREATED)			NOT AVAILABLE		

OPTIONAL			
Slide			
	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm
Oxygen pack			
	OX07880 – 1 UNI SOCKET OX07885 – 2 DIN SOCKETS OX07897 – 2 BS SOCKETS	CODES TO BE CREATED	CODES TO BE CREATED
Electric pack			
	STANDARD EQUIPMENT	STANDARD EQUIPMENT	STANDARD EQUIPMENT
AUTOMATIC HEIGHT ADJUSTMENT SENSOR			
	BO10007	BO10007	BO10007
BASIN FOR HOUSING			
	BO10107	BO10107	NOT AVAILABLE
EXTERNAL CONTROL PANEL			
	BO10105	CODE TO BE CREATED	CODE TO BE CREATED



Model	BO10005 – COMPAS E
General specs	Sideways movable and extractable With tilt front and tilt back position
Surface height	37,5 cm
Loading point when extracted	From 13 to 37,5 cm above installation floor
Sideways movement	230 mm
Extraction tilted	Yes
Suspension	No
Housing	Yes
Double housing version	BO10205

OPTIONAL	
Slide 	BO10100 – 120mm BO10101 – 180mm BO10102 – 350mm BO10104 – 500mm
Oxygen pack 	CODES TO BE CREATED
Electric pack 	STANDARD EQUIPMENT
AUTOMATIC HEIGHT ADJUSTMENT SENSOR 	BO10007
BASIN FOR HOUSING 	BO10107
EXTERNAL CONTROL PANEL 	BO10106