

## CARRERA TEC

### Trolley for incubator



Spencer Carrera Tec is a compact and lightweight self-loading trolley for transport of incubators. The legs fold independently operating the specific commands and open automatically when unloading.

#### Specific Characteristics

The distance between incubator supports is adjustable in order to adapt to your own incubator  
 Integrated stainless steel instruments table is large and easy to clean  
 Color coded commands  
 2 fixed wheels and 2 pivoting with large diameter to make movement easier  
 Brakes on rear wheels

#### Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	32 ± 1 kg
Table dimensions	520x450 mm
Table height <small>(From main frame)</small>	190 mm
Incubator supports distance	From 250 to 800mm
Loading capacity	170 kg
Wheels	ø200 mm
Height of the loading carriage wheel	63 cm
Frame material	Steel

#### Accessories

ST42100A	FP Posterior fastener for self loading stretchers
ST42200A	FA Anterior fastener for self loading stretchers

Class I Med.Device according to Dir. 93/42/CEE

TROLLEY INCUBATOR TEC 63 CHROME  
**CA00010C**



CND Classification

V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala BaganzaPR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)