

# CARRERA XL

## Self loading stretcher



Spencer Carrera XL is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically.

### Specific features

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Patient board with Trendelenburg and Fowler antishock position
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars Safe bar
- Adjustable heights system

### Optional accessories

|          |   |
|----------|---|
| ST42100A | FP Posterior fastener for self loading stretchers |
| ST42200A | FA Anterior fastener for self loading stretchers  |
| EN90003C | END-T Table for stretchers 10g cert.              |
| ST00491A | STX 90 Telescopic headrest for stretchers         |
| ST00497B | DNA STRAP Belt w/integrated re-winding system     |
| ST00498B | DNA STRAP Thorax. Belt w/int.re-winding system    |
| ST00499B | STX 499 - universal adjustable thorax belt        |
| IF01047C | Track 4-30 telescopic IV holder                   |
| IF01049B | Track 5 telescopic IV holder for Carrera          |
| CB09025C | Tanker oxygen tank holder for stretcher side bar  |
| CB09028C | Tanker support for Safe Bar                       |
| ST70000A | QMX777 Black PVC anatomic mattress                |

### Standard equipment

- ST70002A | STX 702 two pieces black reflex metal belts

### Technical data

|                               |         |
|-------------------------------|---------|
| Length                        | 1980 mm |
| Width                         | 585 mm  |
| Weight                        | 37 kg   |
| Loading capacity              | 170 kg  |
| Wheels                        | ø200 mm |
| Height loading carriage wheel | 75 cm   |
| Frame material                | Steel   |
| Side bars length              | 680 mm  |
| Side bars Height              | 200 mm  |

Dimensions and weights are subject to tolerances of ±10mm and ± 1kg

### Class I Med.Device according to Dir. 93/42/CEE

CARRERA XL TF TEC BLACK/B  
YELLOW/F ADJUST. HEIGHTS

CA00015

CND Classification



V0899

Rev.0 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)