

CARRERA TEC

Trolley for incubator



Spencer Carrera Tec is a compact and lightweight self-loading trolley for transport of incubators. The legs fold independently operating the specific commands and open automatically when unloading.

Specific Characteristics

- The distance between incubator supports is adjustable in order to adapt to your own incubator
- Integrated stainless steel instruments table is large and easy to clean
- Color coded commands
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels

Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	32 ± 1 kg
Table dimensions	520x450 mm
Table height <small>(From main frame)</small>	190 mm
Incubator supports distance	From 250 to 800mm
Loading capacity	170 kg
Wheels	ø200 mm
Height of the loading carriage wheel	60 cm
Frame material	Steel

Accessories

ST42100A	FP Posterior fastener for self loading stretchers
ST42200A	FA Anterior fastener for self loading stretchers

Class I Med.Device according to Dir. 93/42/CEE

CARRERA TEC TROLLEY FOR
INCUBATOR
CA10032C



CND Classification	V0899
Registration number	158108

Rev.2 (08/01/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it