

# PRO TROLLEY

## Trolley for incubator



Self-loading trolley for transport of incubators. The legs fold independently operating the specific commands and open automatically when unloading.

### Specific Characteristics

- The distance between incubator supports is adjustable in order to adapt to your own incubator
- Integrated stainless steel instruments table is large and easy to clean
- Color coded commands
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Possibility to place a spinal board in the dedicated compartment

### Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	38 ± 1 kg
Table dimensions	520x450 mm
Table height <small>(From main frame)</small>	190 mm
Loading capacity	150 kg
Wheels	ø200 mm
Height of the loading carriage wheel	60 cm
Frame material	Steel

### Accessories

ST42100A	FP Posterior fastener for self loading stretchers
ST42200A	FA Anterior fastener for self loading stretchers

Class I Med.Device according to Dir. 93/42/CEE  
 PRO TROLLEY HIGH FOR INCUBATOR  
 CA10029C



CND Classification V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)