

CARRERA T TEC

Self loading stretcher



Spencer Carrera T Tec is a compact and lightweight self-loading stretcher for transport of patients. It can be used for rescue operations on road and sanitary vehicles. The legs fold independently operating the specific commands and open automatically when unloading.

Specific features

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars *Safe bar*
- Board with Trendelenburg antishock position

Accessories

| | |
|----------|---|
| ST42100A | FP Posterior fastener for self loading stretchers |
| ST42200A | FA Anterior fastener for self loading stretchers |
| EN90003C | END-T Table for stretchers 10g cert. |
| ST00491A | STX 90 Telescopic headrest for stretchers |
| ST00497B | DNA STRAP Belt w/integrated re-winding system |
| ST00498B | DNA STRAP Thorax. Belt w/integrated re-winding system |
| IF01047C | Track 4-30 telescopic IV holder |
| IF01049B | Track 5 telescopic IV holder for Carrera |
| CB09025C | Tanker oxygen tank holder for stretcher side bar |
| CB09028C | Tanker support for Safe Bar |

Standard equipment

- ST70002A STX 702 two pieces black reflex metal belts

Technical data

| | |
|--------------------------------------|-----------------|
| Length | 1980 mm |
| Width | 580 mm |
| Weight | 37 kg |
| Loading capacity | 170 kg |
| Wheels | ø200 mm |
| Height of the loading carriage wheel | 50 cm |
| Frame material | Stainless steel |
| Side bars length | 680 mm |
| Side bars Height | 200 mm |

Class I Med.Device according to Dir. 93/42/CEE

CARRERA T TEC LOW LOADING
HEIGHT CHROME FRAME BLACK
BOARD



CA10155C

CND Classification

V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it