

CARRERA T ACTIV HIGH

Self Loading stretcher



Spencer Carrera Activ is a self-loading stretcher for transport of patients. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The front part of the frame can be folded, reducing considerably the length of the stretcher when it is needed. The legs fold away independently, activating the dedicated commands during loading procedures, and open automatically when unloading.

Specific Characteristics

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
Board with Trendelenburg antishock position
2 fixed wheels and 2 pivoting with large diameter to make movement easier
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>
Can be shortened in less than 2 seconds
Greater maneuverability in tight spaces than a standard stretcher

Accessories

ST42100A	Posterior fastener for self loading stretchers
ST42200A	Anterior fastener for self loading stretchers
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated re-winding system
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
ST70000A	QMX777 Black PVC anatomic mattress
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard Equipment

ST70002A	STX 702 two pieces black reflex metal belts
ST70018A	QMX 776 BLACK MATTRESS, PVC, 1 SECTION

Technical data

Length	1975 ± 10 mm
Width	580 ± 10 mm
Length shortened	1680 ± 20 mm
Weight	36 ± 1 kg
Loading capacity	170 kg
Wheels	ø200 mm
Height loading carriage wheel	60 cm
Frame material	Stainless steel
Patient board height	
Lowered stretcher	280 ± 30 mm
Raised stretcher	870 ± 30 mm
Height at the top of the side bar	
Lowered stretcher	480 ± 30 mm
Raised stretcher	1070 ± 30 mm
Side bars length	680 mm
Possible heights	2

Class I Med.Device according to Dir. 93/42/CEE

CARRERA T ACTIV HIGH CHROME
FRAME BLACK BOARD

CA20052C

CND Classification



V0899

Rev.1 (18/07/17)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it