

SUNNY STRETCHER

Self-loading stretcher



The product can be subject to changes without notice



Sunny stretcher is a compact and lightweight self-loading stretcher for transport of patients. It can be used for rescue operations on road and sanitary vehicles.

Specific Characteristics

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars

Technical data

Length	1805 ± 10 mm
Width	570 ± 10 mm
Weight	38 ± 1kg
Loading capacity	250 kg
Wheels	∅ 150 mm
Height of the loading carriage wheel	63 cm
Frame material	Al, Steel

Class I Med.Device according to Dir. 93/42/CEE

MOD.S041 CARRERA SUNNY
STRETCHER

CA31054S

CND Classification



V0899

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it