

CARRERA XL

Barella autocaricante 10g



Spencer Carrera XL is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. 10G certified when used with the fixing system Spencer R-Max.

Specific features

| |
|---|
| Polyethylene board easy to clean |
| Color coded commands |
| Adjustable backrest |
| Antishock positions Trendelenburg |
| 2 fixed wheels and 2 pivoting with large diameter to make movement easier |
| Brakes on rear wheels |
| Foldable side bars Safe bar |
| Telaio verniciato giallo |
| Adjustable heights system to facilitate patient transfer procedures |

Accessories

| | |
|----------|---|
| ST42702C | R-Max fixing system <small>(Mandatory for 10G versions)</small> |
| EN90003C | END-T Table for stretchers 10g cert. |
| ST00491A | STX 90 Telescopic headrest for stretchers |
| ST00497B | DNA STRAP Belt w/integrated re-winding system |
| ST00498B | DNA STRAP Th.Belt w/integr.re-winding system |
| IF01047C | Track 4-30 telescopic IV holder |
| IF01049B | Track 5 telescopic IV holder for Carrera |
| ST70000A | QMX777 Materasso termosaldato anatomico PVC nero |
| CB09025C | Tanker oxygen tank holder for stretcher side bar |
| CB09028C | Tanker support for Safe Bar |

Standard Equipment

| | |
|----------|---|
| ST70002A | STX 702 two pieces black reflex metal belts |
|----------|---|

Technical data

| | |
|-------------------------------|-----------------|
| Length | 1980 ± 10mm |
| Width | 585 ± 10mm |
| Weight | 37 ± 1 kg |
| Loading capacity | 170 kg |
| Wheels | ø200 mm |
| Height loading carriage wheel | 75 cm |
| Frame material | Stainless steel |
| Side bars length | 680 mm |
| Side bars Height | 200 mm |

Class I Med.Device according to Dir. 93/42/CEE

CARRERA XL T TEC BLACK BOARD
YELLOW FRAME WITH ADJUSTABLE
HEIGHTS 10G



CA60000C

| | |
|---------------------|--------|
| CND Classification | V0899 |
| Registration number | 789248 |

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it