

CARRERA CINCO

10 g Self loading stretcher with 5 heights



Spencer Carrera Cinco is a self-loading stretcher with 5 adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. 10G certified when used with the fixing system Spencer R-Max or S-MAX.

Specific features

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest in 5 position
- Antishock positions Trendelenburg
- Extreme maneuverability thanks to the Twist system which, when operated, makes pivoting also the front wheels.
- Brakes both on rear and front wheels
- Foldable side bars *Safe bar*
- 5 Adjustable heights system

Optional accessories

Code	Description
ST42707B	S-Max fastener (Mandatory for 10g versions)
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Th.Belt w/integr.re-winding system
ST00499B	STX 499 – 4 pcs adjustable thorax belt
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard Equipment

- 3 belts with quick release metal buckle
- Mattress

Technical data

Length	1970 mm
Width	570 mm
Patient floor height when folded	Approx. 30cm
Weight	43 Kg
Loading capacity	200 kg
Wheels	∅200mm ± 5%
Height loading carriage wheel	68 cm
Frame material	Stainless steel
Side bars length	680 mm
Side bars Height	200 mm

Dimensions and weight subject to tolerances of ±10mm and ± 3kg

Class I Med.Device according to Dir. 93/42/CEE

MOD.Z001 CARRERA CINCO T
HIGH 68 F/CHROM B/BLACK 10G
CA60005



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it