

# CARRERA TEC

## Trolley for incubator



Spencer Carrera Tec is a compact and lightweight self-loading trolley for transport of incubators. The legs fold independently operating the specific commands and open automatically when unloading.

### Specific Characteristics

- The distance between incubator supports is adjustable in order to adapt to your own incubator
- Integrated stainless steel instruments table is large and easy to clean
- Color coded commands
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels

### Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	33 ± 1 kg
Table dimensions	520x450 mm
Table height <small>(From main frame)</small>	190 mm
Incubator supports distance	From 250 to 800mm
Loading capacity	150 kg
Wheels	ø200 mm
Height of the loading carriage wheel	60 cm
Frame material	Steel

### Accessories

ST42702C	R-Max fixing system <small>(Mandatory)</small>
----------	--

### Class I Med.Device according to Dir. 93/42/CEE

YELLOW TROLLEY TEC HIGH FOR  
INCUBATOR 10G  
CA60022C



### Compliant with

EN 1789 for the coupling between the fastener and the trolley  
**Note:** Compliance with the requirements of EN 1789 for fixations between the incubator and the trolley is an installer responsibility

CND Classification	V0899
Registration number	343854

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)