

## CARRELLO TEC

### Trolley for incubator



Spencer Carrera Tec is a compact and lightweight self-loading trolley for transport of incubators. The legs fold independently operating the specific commands and open automatically when unloading.

#### Specific features

The distance between incubator supports is adjustable in order to adapt to your own incubator

Integrated stainless steel instruments table is large and easy to clean

Color coded commands

2 fixed wheels and 2 pivoting with large diameter to make movement easier

Brakes on rear wheels

#### Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	32 ± 1kg
Table dimensions	520x450 mm
Table height <small>(From the main frame)</small>	190 mm
Loading capacity	150 kg
Wheels	ø200 mm
Height of the loading carriage wheel	60 cm
Frame material	Steel

#### Accessories

ST42707B 10G fastener S-MAX (Mandatory)

#### Compliant with

EN 1789 for the coupling between the fastener and the trolley

**Note:** Compliance with the requirements of EN 1789 for fixations between the incubator and the trolley is an installer responsibility

Class I Med.Device according to Dir. 93/42/CEE

TEC TROLLEY CHROME HIGH  
VERSION FOR INCUBATOR 10G

CA60023C

CND Classification



V0899

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)