

CARRERA XL

10g Self Loading stretcher



Spencer Carrera XL is a self-loading stretcher with adjustable heights for patient transport. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. 10G certified when used with the fixing system Spencer S-Max.

Specific features

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Antishock positions Trendelenburg and Fowler
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars *Safe bar*
- Polished frame
- Adjustable heights system
- System to lock the legs in closed position

Accessori opzionali

ST42707B	S-Max fixing system (Mandatory for 10g versions)
ST70000A	QMX777 PVC anatomic mattress
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Th.Belt w/integr.re-winding system
ST00499B	STX 499 – 4 pcs adjustable thorax belt
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard Equipment

- ST70002A | STX 702 two pieces black reflex metal belts

Technical data

Length	1980 mm
Width	585 mm
Weight	36 kg
Loading capacity	170 kg
Wheels	ø200 mm
Height loading carriage wheel	75 cm
Frame material	Stainless Steel
Side bars length	680 mm
Side bars Height	200 mm

Dimensions and weight subject to tolerances of ±10mm and ± 1kg

Class I Med.Device according to Regulation 2017/745/UE

CARRERA XL TF TEC BLACK BOARD
CHROME FRAME WITH
ADJUSTABLE HEIGHTS 10G

CA60026C

ClassificationCND



V0899

Rev.0 (24/03/2022)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it