

**CARRERA TEC**

Self loading stretcher



A wholesome and solid stretcher, able to rewrite the parameters of excellence for its category. No one before has ever given such care and attention to technological innovation and exploitation of space.. Designed, developed and constructed entirely in Europe to satisfy all the needs of those who in the rescue field, are used to high standards.

The Carrera stretcher is supplied with two straps model STX 702.

10G Certified if used with R-MAX fixation system.

**TECHNICAL DATA**

Length	1970 mm
Width	570 mm
Loading height mm	From 550 mm to 690
Wheels	Ø200 mm
Brakes	On Rear wheels
Frame material	Steel
Patient floor material	PE black
Patient floor type	Trendelenburg
Load capacity	170 kg
Weight	32 kg
Side bars height	200 mm
Side bars length	680 mm

CA60154B

CARRERAT TEC HIGH  
F/BLACK F/CHROME  
10G CERT.



UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>