

## CARRERA TEC

### 10g Self Loading Stretcher



Spencer Carrera T Tec is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during the unloading procedure, they open automatically. 10G certified if used with dedicated fastener S-MAX.

#### Specific features

Polyethylene board easy to clean  
 Color coded commands  
 Adjustable backrest (from 0° to 75°)  
 Trendelenburg antishock position  
 2 fixed wheels and 2 pivoting with large diameter to make movement easier  
 Brakes on rear wheels  
 Foldable side bars *Safe bar*

#### Optionals

ST42707B	S-Max fastener (Mandatory for 10G versions)
ST70000A	PVC anatomic mattress QMX/777
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax.Belt w/integr.re-winding syst.
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

#### Standard equipment

STX 702 two pieces black reflex metal belts

#### Technical data

Length	1970 mm
Width	570 mm
Weight	38 kg
Loading capacity	170 kg
Wheels	ø200 mm
Height loading carriage wheel	60 cm
Frame material	Stainless steel
Side bars length	680 mm
Side bars Height	200 mm

Dimensions and weight subject to tolerances of ±10mm and ± 1kg

Class I Med. Dev. according to Regulation 2017/745/UE

CARRERA TEC T HIGH LOADING  
 BLACK/B. CHROME/F. 10 G  
**CA60154B**



**CND Classification**

V0899

Rev.1 (05/11/2021)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)