

## CARRERA TF TEC HIGH

### 10G Self Loading stretcher



Spencer Carrera TF Tec is a self-loading stretcher for transport of patients. It is a compact and light stretcher, usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands during loading procedures, and open automatically when unloading. 10G Certified if used with Spencer S-MAX certified fastener.

#### Specific Characteristics

Polyethylene board easy to clean  
 Color coded commands  
 Adjustable backrest  
 2 fixed wheels and 2 pivoting with large diameter to make movement easier  
 Brakes on rear wheels  
 Foldable side bars Safe bar  
 Board with Trendelenburg and Fowler antishock positions

#### Accessories

ST42707B	S-Max fastener <small>(Mandatory for 10G versions)</small>
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integrated reel
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

#### Standard Equipment

ST70002A | STX 702 two pieces black reflex metal belts

#### Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	38 ± 1 kg
Loading capacity	170 kg
Wheels	ø200 mm
Height loading carriage wheel	60 cm
Frame material	Stainless steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Dir. 93/42/CEE

CARRERA TF TEC BLACK BOARD  
 YELLOW FRAME 10G. CERT.

CA61254C

**CND Classification** V0899  
**Registration number** 160512



Rev.0 (22/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)