

CARRERA S TEC 75

Self-loading stretcher



Spencer Carrera S Tec is a compact and lightweight self-loading stretcher for transport of patients. It can be used for rescue operations on road and sanitary vehicles. The legs fold independently operating the specific commands and open automatically when unloading. 10G certified if used with the dedicated fastener

Specific Characteristics

Polyethylene board easy to clean
Color coded commands
Adjustable backrest
2 fixed wheels and 2 pivoting with large diameter to make movement easier
Brakes on rear wheels
Foldable side bars <i>Safe bar</i>
Yellow painted frame

Accessories

ST42707B	S-Max fastener <small>(Mandatory for 10G versions)</small>
ST42702C	R-Max fastener <small>(Mandatory for 10G versions)</small>
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integr.re-wind.syst.
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

Standard Equipment

ST70002A	STX 702 two pieces black reflex metal belts
ST70018A	PVC anatomic mattress 1 piece

Technical data

Length	1980 ± 10mm
Width	580 ± 10mm
Weight	37 ± 1kg
Loading capacity	170 kg
Wheels	ø200 ± 5% mm
Height of the loading carriage wheel	75 cm
Frame material	Steel
Side bars length	680 mm
Side bars height	200 mm

Class I Med.Device according to Dir. 93/42/CEE

CARRERA S TEC 75 F/YELLOW
P/BLACK W/QMX776 10G CER
CA70021



CND Classification

V0899

Rev.0 (03/12/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it