

# CARRERA T TEC MAX

## Self-loading stretcher



Spencer Carrera T Tec Max is a compact and lightweight self-loading stretcher for transport of patients. It can be used for rescue operations on road and sanitary vehicles. The legs fold independently operating the specific commands and open automatically when unloading. 10G certified if used with the dedicated fastener.

### Specific Characteristics

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars *Safe bar*
- Yellow painted frame
- Trendelenburg Antishock position

### Accessories

ST42707B	S-Max fastener (Mandatory)
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integ. re-winding syst.
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
ST70000A	QMX777 Black PVC anatomic mattress
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar

### Standard Equipment

- STX 702 two pieces black reflex metal belts
- QMX777 PVC anatomic mattress

### Technical data

Length	1980 ± 10 mm
Width	580 ± 10 mm
Weight	38 ± 1kg
Loading capacity	250 kg
Wheels	ø200 mm
Height of the loading carriage wheel	68 cm
Frame material	Steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Regulation 2017/745/UE

CARRERA T TEC 68 MAX  
F/CHROME P/BLACK 10G

CA70022

CND Classification



V0899

Rev.0 (29/11/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)