

CARRERA TEC LOW

10 Self Loading stretcher



Spencer Carrera T Tec is a compact and lightweight self-loading stretcher for transport of patients. It can be used for rescue operations on road and sanitary vehicles. The legs fold independently operating the specific commands and open automatically when unloading. 10G certified if used with the dedicated fastener.

Specific features

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest
- Trendelenburg antishock position
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars *Safe bar*

Accessories

ST42707B	S-Max fastener (Mandatory for 10G versions)
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integr.re-wind.syst.
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar
EN90003C	END-T Table for stretchers 10g cert.

Standard equipment

ST70018A	QMX776 PVC anatomic mattress
ST70002A	STX 702 two pieces black reflex metal belts

Technical data

Length	1980 mm
Width	585 mm
Weight	35 kg
Loading capacity	170 kg
Wheels	ø200 mm
Height of the loading carriage wheel	50 cm
Frame material	Stainless steel
Side bars length	680 mm
Side bars height	200 mm

Dimensions and weight subject to tolerances of ±10mm and ± 1kg

Class I Med. Dev. according to Reg. 2017/745/UE

LOW CARRERA T TEC P/BLACK
F/CHR. W/QMX776 10G CER.

CA70051C

CND Classification



V0899

Rev.0 (01/09/2022)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it