

CARRERA PRO MAX

Self-Loading stretcher

With housing for spine board



Carrera Pro is a self-loading stretcher for patient transport usable in rescue operations on streets and rescue vehicles. The legs fold away independently, activating the dedicated commands and they open automatically during the unloading procedure. The spinal board can be placed in a dedicated housing located under the patient board. 10G Certified if used with Spencer S-MAX certified fastener.

Specific features

- Polyethylene board easy to clean
- Color coded commands
- Adjustable backrest (from 0° to 75°)
- Antishock position Trendelenburg
- 2 fixed wheels and 2 pivoting with large diameter to make movement easier
- Brakes on rear wheels
- Foldable side bars *Safe bar*
- Spinal board easy to remove, making immobilization procedures faster

Accessories

ST42707B	S-Max fastener <small>(Mandatory for 10G versions)</small>
EN90003C	END-T Table for stretchers 10g cert.
ST00491A	STX 90 Telescopic headrest for stretchers
ST00497B	DNA STRAP Belt w/integrated re-winding system
ST00498B	DNA STRAP Thorax. Belt w/integr. re-wind. syst.
IF01047C	Track 4-30 telescopic IV holder
IF01049B	Track 5 telescopic IV holder for Carrera
CB09025C	Tanker oxygen tank holder for stretcher side bar
CB09028C	Tanker support for Safe Bar
ST02061B	B-BAK PIN YELLOW SPINE BOARD W/PINS

Note

Is possible to use only the spinal board Spencer B-Bak PIN not included

Standard Equipment

ST70002A	STX 702 two pieces black reflex metal belts
ST70018A	QMX776 Anatomic mattress with hermetical seamless design, antibacterial and fireproof, resistant to disinfectants

Technical data

Length	1970 mm
Width	570 mm
Weight	38 ± 1 kg
Loading capacity	250 kg
Wheels	ø200 mm
Height loading carriage wheel	60 cm
Frame material	Stainless steel
Side bars length	680 mm
Side bars Height	200 mm

Class I Med.Device according to Dir. 93/42/CEE

HIGH CARRERA T PRO F/CH. P/BLACK
W/QMX776 10G CER.

CA71152C



CND Classification

V0899

Rev.0 (10/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it