

## ALFA 2

### Multipurpose backpack



Multipurpose emergency backpack extremely versatile and robust. It is characterized by a large volume body, with numerous compartments and equipped with various internal fixation elements

#### Specific features

Carrying handle and adjustable shoulder straps to make carrying easier and increase comfort

PVC coated bottom with TPU support elements, offer great protection and make cleaning procedures easier

2 side pockets with zip closure and elastic bands for fixation of small objects

Main body divided into 3 compartments with "book" type opening

Each of the 3 compartments is characterized by a rigid perimetral reinforcement which provides protection to the content and improves stability

Wide reflective inserts to increase visibility

Transparent front pocket with strap closure

Front pocket in the upper part with strap closure

6 fixed pockets in the upper compartment to place accessories

5 fixation belts with strap closure on 2 sides of the upper compartment

On the internal side of the upper closure flap, there is foldable fabric strip with elastic bands and at which end there is a further pocket with strap closure

#### Technical data

Dimensions	630x450x250 ± 30 mm
Weight	2,35 ± 0,2 kg
Volume	43 l
Materials	Nylon, PE, TPU, Polyester

ALFA 2 - MULTI-PURPOSE EMERGENCY BACKPACK  
CB00300B

Rev.0 (22/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)