


MARK PRO

Complete professional ALS backpack



Professional backpack characterized by extreme versatility thanks to the arrangement of the internal parts, to the protection offered to the content and to the dual transport possibility.

Specific features

- Equipped with two handles to allow transport also as a bag. The handles are both on the upper part and on the side
- 1 outside pocket with zip closure on the bottom and 1 frontal pocket with access from the side
- 2 outside pockets with zip closure for mobile radio or phone
- 3 transparent pockets with zip closure on the internal side of the main compartment
- Internal partition with 2 big pockets with zip closure. One is transparent
- 2 transparent pockets with zip closure on the internal side of the frontal flap
- Backrest reinforced with aluminum structure
- Padded main compartment to provide protection to the content
- PVC covered bottom to minimize accumulation of dirt and to make cleaning easier
- Padded shoulder straps adjustable in height, with chest and abdominal belt with quick release buckle for a better load distribution on the operator's back
- Backrest with paddings and covered with breathable fabric
- Reflective inserts to increase visibility

Equipment

BG02027	PRO-B-LIFE PVC ADULT RESUSCITATOR
BS04011	2x BURN S-10X10 STERILE BURN STOP DRESSING
BS04021	2x BURN S- 20X20 STERILE BURN STOP DRESSING
CB03030	PHARMA TANK 1 SEMIRIGID PHIAL HOLDER
DG00150	YANKEUR – STERILE RIGID SPIRATION CANNULA
DG00212	4x DGX 50 CH/FG 12 SUCTION CATHETER
DG00603	UNIVERSAL BLACK TUFF CUT SCISSORS
DG01410	ADULT SPHYGMOMANOMETER
DG02110	2x DGX 110 BLACK PHONENDOSCOPE
DG08062	DGX 062 PEN LIGHT WITHOUT BATTERIES
DG08005	IME COMPLETE INTUBATION KIT
DG03100	LARYNGOSCOPE HANDLE TRADITIONAL LIGHT
DG05001	MACINTOSH BLADE SZ 1
DG05002	MACINTOSH BLADE SZ 2
DG05003	MACINTOSH BLADE SZ 3
DG80550	CUFFED STERILE TRACHEAL TUBE SZ 5
DG80560	CUFFED STERILE TRACHEAL TUBE SZ 6
DG80570	CUFFED STERILE TRACHEAL TUBE SZ 7
DG80580	CUFFED STERILE TRACHEAL TUBE SZ 8
FG00605	2x MAGILL ADULT FORCEPS 268 mm
FG00607	2x MAGILL PEDIATRIC FORCEPS 200 mm
GM03050	GMX 05 STERILE DISPOSABLE SYRINGE 5ML
GM03215	S PLASTER TAPE ADHESIVE PLASTER 1,25cm x 9,14m
RM10320	GUEDEL DISPOSABLE CANNULA 50 mm
RM10321	GUEDEL DISPOSABLE CANNULA 60 mm

RM10322	GUEDEL DISPOSABLE CANNULA 70 mm
RM10323	GUEDEL DISPOSABLE CANNULA 80 mm
RM10324	GUEDEL DISPOSABLE CANNULA 90 mm
RM30700	MOUNDER JAW SPREADER
FG03308	GMX 308 SURGICAL FORCEPS 140 mm
FG08810	EXELL SILICONE MOUTH WEDGE
GM00363	2x GMX 363 KLEMMER FORCEPS 163mm
GM02215	T-PLASTER REEL OF PLASTER 1,25 cm x 5 m
GM02225	T-PLASTER REEL OF PLASTER 2,5 cm x 5 m
GM02250	T-PLASTER REEL OF PLASTER 5 cm x 5 m
GM02405	2x G- GAUZE BENDAGE 5 cm x 5 m
GM02407	2x G- GAUZE BENDAGE 7 cm x 5 m
GM02410	2x G- GAUZE BENDAGE 10cm x 5 m
GM03050	2x GMX05 STERILE DISPOSABLE SYRINGE 5 ml
GM03100	2x GMX10 STERILE DISPOSABLE SYRINGE 10 ml
GM04020	SSP – PLASTERS 40PCS
GM04250	GMX 003 DISINFECTANT 250 ml
GM04300	GMX 300 HYDROGEN PEROXIDE 10VOL 250 ml
IF00100	2x TIV AUTOMATIC TOURNIQUETTE
IF00421	2x SPENCER ICE PAK TNT
RM10310	2x KIT GUEDEL 6 SZ W/BOX
RM10889	3x CUFFED MASK SZ 0 W/OUT VALVE
RM10893	CUFFED MASK SZ 4 W/OUT VALVE
RM30300	COLLIN STAINLESS STEEL TONGUE FORCEPS
SC75100	JET COMPACT 300D
SC75011	TRANSPORT BAG FOR JET COMPACT
SC75009	SHOULDER STRAP FOR JET
ST00630	2x HT ADULT ISOTHERMIC SHEET
TA09020	FLEX MASK W/BOX
TT04020	3x SPENCER STICK YELLOW LUMINOUS STICK

Technical data

Dimensions	360x520x300±10 mm
Empty backpack weight	3,13 ± 0,03 kg
Backpack volume	38 ± 1 l
Loading capacity	30 kg
Materials	Nylon, PE, PVC

MARK COMPLETE ALS BACKPACK BLACK/YELLOW
CB00428A

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it