

A911

Rescue backpacks set



A911 is a set of backpacks extremely practical and versatile thanks to their compactness and lightness. The single backpack features a padded main compartment to ensure the safety of the content and has a reinforced structure to limit deformation. The bright coloring and reflective strips increase its visibility. Numerous internal and external pockets, allow an easier organization of the space. The bottom is reinforced and non-slip as well as the carrying handle.

The sales code includes 6 backpacks.

TECHNICAL DATA

Dimensions	490x470x210 ± 30 mm
Weight	1,97 ± 0,1
Materials	Nylon and TPU
External pockets	4
Internal bags	5 detachable with zip closure and transparent top window
Dimensions of an internal bag	280x120x70 ± 10 mm
Portability	Carrying handle and padded shoulder straps

CB00911A

A911 SET 6 RED BACKPACKS
W/5 BAGS

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>