

A911

Complete BLS Backpack



Compact professional backpack equipped with devices most commonly used in BLS procedures.

Included Equipment		Specific Characteristics	
BG02007A	PRO B-LIFE - PVC AD. RES. W/OVERP.VALVE R-VALVE	Dual transport possibility thanks to the padded shoulder straps and the soft carrying handle placed on the upper side	
BS04011D	BURN S (BA) 10 X 10 CM - BURN STOP DRESSING	Adjustable shoulder straps with adjustable thoracic and abdominal cross bands with quick release buckle to provide more stability during transport	
BS04021D	BURN S (BA) 20 X 20 CM - BURN STOP DRESSING	High visibility thanks to the reflective bands	
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS	The frontal pocket is capacious and can be completely opened. On both internal sides are sewn elastic bands to attach accessories	
DG01410C	DGX 411 - AN.SPHIGMO. 2TUBES W/AD. CUFF BLUE	Large pocket outside the frontal flap equipped with zip closure	
DG02110B	ECONOMICAL PHONENDOSCOPE	The main compartment is padded and is characterized by	
DG08033D	TWK33 ALUM. TORCH W/OUT BATT	- elastic bands on the perimeter suitable for accessories storage	
FG00605A	FORCEPS MAGILL ADULT	- zip closure that allow the complete opening of the compartment	
FG00607A	MAGILL PAEDIATRIC - TUBE-GUIDING FORCEPS	- strap bands to keep in place the internal pockets	
FG03308B	GMX 308 - SURGICAL FORCEPS	- perimetral reinforcement to limit deformations	
FG08810A	MOUTH WEDGE EXCELL	2 wide lateral pockets with zip closure	
GM00363G	KLEMMER FORCEPS	5 internal detachable pockets with zip closure and transparent top	
GM00373G	KLEMMER CURVED FORCEPS	Non-slip rubber coated bottom to make cleaning easier	
GM02215K	REEL OF PLASTER 1,25 cm x 5 m	Technical data	
GM02225K	ADHESIVE BANDAGE ROLL FABRIC 2,5CMX5M	Dimensions	490x470x210 ± 30 mm
GM02250K	CLOTH PLASTER 5CM X 5M	Dimensions of each internal pocket	280x120x70 ± 10 mm
GM02405K	GAUZE BENDAGE 5M x5CM	Materials	Nylon, TPU, PVC
GM02407K	GAUZE BANDAGE 5M X 7CM	Volume	34 ± 1 l
GM02410K	GAUZE BANDAGE 5M X 10CM	Weight of empty backpack w/pockets	1,97 ± 0,1 kg
GM03100K	DISPOSABLE STERILE SYRINGE 10 ML	Medical devices CE marked according to dir. 93/42/CE	
GM04020K	PLASTERS - 40 PIECES	A911 - COMPLETE BLS BACKPACK W/5 POCKETS	
GM04250K	GMX 003 – DISINFECTANT 250 ml	CB00913A	
GM04300K	HYDROGEN PEROXIDE 10 VOL. 250 ml		
IF00100C	T.I.V. SPENCER AUTOMATIC TOURNIQUETTE		
IF00421A	COLD PACK SPENCER ICE TNT SINGLE USE BOX 25 Pcs		
RM10311A	GUEDEL AIRWAY KIT 8 SIZES W/BOX		
RM10881A	CUFFED MASK PVC SZ.1 AIR CUSHION		
RM10882A	CUFFED MASK PVC SZ.2 AIR CUSHION		
RM10883A	CUFFED MASK PVC SZ.3 AIR CUSHION		
RM10884A	CUFFED MASK PVC SZ.4 AIR CUSHION		
RM10885A	CUFFED MASK PVC SZ.5 AIR CUSHION		
RM30300A	COLLIN - STAINLESS STEEL TONGUE FORCEPS		
ST00630A	ADULT ISOTHERMIC SHEET		
TA09020A	FLEX MASK – PVC MASK W/BOX		
TT04010E	RED CYALUME LUMINOUS STICK		
TT04020E	YELLOW CYALUME LUMINOUS STICK		
TT04030E	ORANGE CYALUME LUMINOUS STICK		
TT04040E	GREEN CYALUME LUMINOUS STICK		

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it