

A911

Emergency backpack with bags



Professional compact and light backpack designed to accommodate devices most commonly used in BLS procedures.

Specific features

Dual transport possibility thanks to the padded shoulder straps and the soft carrying handle placed on the upper side

Adjustable shoulder straps with adjustable thoracic and abdominal cross bands with quick release buckle to provide more stability during transport

High visibility thanks to the reflective bands

The frontal pocket is capacious and can be completely opened. On both internal sides are sewn elastic bands to attach accessories

Large pocket outside the frontal flap equipped with zip closure

The main compartment is padded and is characterized by

- elastic bands on the perimeter suitable for accessories storage
- zip closure that allow the complete opening of the compartment
- strap bands to keep in place the internal pockets
- perimetral reinforcement to limit deformations

2 wide lateral pockets with zip closure

5 internal detachable pockets with zip closure and transparent top

Non-slip rubber coated bottom to make cleaning easier

Technical data

Dimensions	490x470x210 ± 30 mm
Dimensions of each internal pocket	280x120x75 ± 10 mm
Volume	34 ± 1L
Materials	Nylon, TPU, PVC
Weight of empty backpack w/pockets	1,95 ± 0,1 kg

A911 - BLUE BACKPACK CORDURA W/5 INTERNAL POUCHES
CB00915A

Rev.0 (01/10/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it