

## T SERIES

### T1



PVC bag suitable to store devices commonly used during intubation procedures.

#### Specific Characteristics

Semi transparent to allow view of the content  
 Zip closure for easy access to the content  
 2 female strap bands to allow attachment inside a bag equipped with male straps

#### Technical data

Dimensions	330x110x30 mm
Weight	50 g
Material	PVC

T1  
 CB02070B

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)