

PHARMA TANK 2

Semi-rigid phial holder



Compact and versatile semi-rigid phial holder with transport handle and zip closure.

Specific features

The case provides excellent protection to the content and is equipped with two internal detachable separator where is possible to store phials of different sizes
 Zip closure for easy access to the content
 2 mesh pockets with elastic band on one side of the separator
 Reflective insert to increase visibility

Technical data

Dimensions	210x210x75 ± 5 mm
Weight	320 ± 20 g
Phial holders	27
Material	EVA

PHARMA TANK 2 - SEMI-RIGID PHIAL
 HOLDERS
CB03031A

Rev.3 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it