

# REUSS

## Back waist bag



Multipurpose pouch realized with extremely robust fabric and with separate compartments for a perfect management of space.

### Specific Characteristics

- The handle on the top allows to carry it as a bag
- Pockets and compartments with zip closure to allow quick access to the content
- Central compartment with frontal opening and elastic band to store up to 7 phials
- Flap of the central compartment with pocket equipped with elastic band to store up to 9 phials
- Lateral compartments with access from top to provide maximum accessibility
- Removable keychain and little pocket on the belt
- Padded belt coated with perforated fabric to increase comfort and breathability
- Belt closure with quick release buckle
- Coated and reinforced bottom

### Technical data

Dimensions <sup>1</sup>	150x210x400±10 mm
Circumference	From 93 to 120 ± 10 cm
Overall volume	7 ± 0,3 l
Weight	490 ± 10 g
Materials	Nylon, PVC

<sup>1</sup> Belt excluded

REUSS - BACK WAIST BAG  
**CB03131A**

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)