

# MAXIZAINO

## Professional EMS backpack



Capacious professional backpack designed to accommodate all the necessary equipment used in ALS procedures.

Included Equipment		Specific Characteristics	
BG02007A	PRO B-LIFE - PVC AD. RES. W/OVERP.VALVE R-VALVE	Dual transport possibility thanks to the padded shoulder straps and carrying handle	
BS04011D	BURN S (BA) 10 X 10 CM - BURN STOP DRESSING	Detachable and adjustable shoulder straps with adjustable cross band to provide more stability during transport	
BS04021D	BURN S (BA) 20 X 20 CM - BURN STOP DRESSING	High visibility thanks to the reflective bands	
CB03030A	PHARMA TANK 1 - SEMI-RIGID PHIAL HOLDERS	2 capacious front pockets with zip closure	
DG00150A	YANKAUER PROBE	Detachable semirigid inner compartment equipped with zip closure and pocket with strap closure	
DG00212C	CATHETER CONTROLLED SUCTION CH/FG 12	In the main body of the backpack are present:	
DG00603A	UNIVERSAL BLACK - TUFF CUT SCISSORS	- two side pockets with zip closure	
DG01411C	DGX 411 - AN.SPHIGMO. 2TUBES W/AD. CUFF ORAN.	- one pocket with zip closure on the inner side of the front flap and a pocket with strap closure	
DG02110B	ECONOMICAL PHONENDOSCOPE	- adjustable straps to fix a cylinder up to 2l	
DG08062D	DGX 062 - BLUE PLASTIC PEN LIGHT NO BATTERIES	- strap bands for attachment of internal pockets	
DG80005A	IME INTUBATION KIT COMPLETE		
FG00605A	FORCEPS MAGILL ADULT		
FG00607A	MAGILL PAEDIATRIC - TUBE-GUIDING FORCEPS		
FG03308B	GMX 308 - SURGICAL FORCEPS		
FG08810A	MOUTH WEDGE EXCELL		
GM00363G	KLEMMER FORCEPS		
GM02215K	REEL OF PLASTER 1,25 cm x 5 m		
GM02225K	ADHESIVE BANDAGE ROLL FABRIC 2,5CMX5M		
GM02250K	CLOTH PLASTER 5CM X 5M		
GM02405K	GAUZE BENDAGE 5M x5CM		
GM02407K	GAUZE BANDAGE 5M X 7CM		
GM02410K	GAUZE BANDAGE 5M X 10CM		
GM03050K	STERILE DISPOSABLE SYRINGE 5ML		
GM03100K	DISPOSABLE STERILE SYRINGE 10 ML		
GM04020K	PLASTERS - 40 PIECES		
GM04250K	GMX 003 – DISINFECTANT 250 ml		
GM04300K	HYDROGEN PEROXIDE 10 VOL. 250 ml		
IF00100C	T.I.V. SPENCER AUTOMATIC TOURNIQUETTE		
IF00421A	COLD PACK SPENCER ICE TNT SINGLE USE BOX 25 Pcs		
RM10310A	GUEDEL AIRWAY KIT 6 SIZES		
RM10880A	CUFFED MASK PVC SZ.0 AIR CUSHION		
RM10893A	CUFFED MASK SIZE. 4 PVC CO. COD WITHOUT VALVE		
RM30300A	COLLIN - STAINLESS STEEL TONGUE FORCEPS		
SC75011E	TRANSPORT BAG FOR JET COMPACT		
SC75100D	JET COMPACT 300 D		
ST00630A	ISOTHERMIC SHEET ADULT SIZES 160X220 CM.		
TA09020A	FLEX MASK SET WITH ORANGE PVC CASE		
		<b>Technical data</b>	
		Dimensions	500x420x220 mm
		Detachable compartment dimensions	395x140x200 mm
		Weight of the empty backpack	3,19 kg
		<b>Suction device technical data</b>	
		Suction flow (free flow)	32l/min
		Vacuum range	From 50 to 550 mmHg
		IP classification	IP12
		Weight	1,6 kg
		<b>Medical devices CE marked according to dir. 93/42/CE</b>	
		MAXIZAINO COMPLETE ALS-PROFESSIONAL EMS BACKPACK	
		CB04023A	

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)